

Every time before printing, make sure that the build plate is clockwise locked.
When the build plate is in place, the blue mark on the knob should show up in front of you.



● Orange ● Grey ● Ivory Net Weight: 1 kg

Ortho Model Resin

Ideal for Orthodontic Models

RAYSHAPE Ortho Model combines accuracy and precision to set brackets with flexibility accurately. It is ideal for orthodontic models. Using Ortho Model dental resin will save you time and labor for a more efficient workflow.

Benefits

- ✔ Precise
- ✔ Flexible
- ✔ Accurate
- ✔ Low Shrinkage

Applications

- U Splints
- U Orthodontic model
- U Aligners

ROI

Amount: 8.9 g/set
93 printed model with a bottle of 1KG



Viscosity

Lower viscosity values indicate that the resin is easier to mix.

387 mPa·s (25°C)

Flexural Strength

Higher flexural strength ensures that the model doesn't break during regular use.

89 MPa

Flexural Modulus

Higher flexural modulus prevents the model from distorting, which can impact accuracy.

3050 MPa